



ECD01
Suction Only



ECD03
Enteral Feeding + Suction

External Heat Exchanger	Stryker Altrix Temperature Management System or Gaymar/Stryker Medi-Therm III	Stryker Altrix Temperature Management System or Gaymar/Stryker Medi-Therm III
Heat Exchange Connector	Gaymar/Stryker Klik-Tite	Gaymar/Stryker Klik-Tite
Enteral Feeding Connector	(none)	ENfit
Intended Duration of Use	72 h in the US / 120 h outside the US	72 h in the US / 120 h outside the US
Material	Medical-grade silicone	Medical-grade silicone
Water Volume	55 mL (1.9 fl oz.)	55 mL (1.9 fl oz.)
Outside Diameter	12.0 mm (0.47 in, 36 Fr)	12.0 mm (0.47 in, 36 Fr)
Inside Diameter of Gastric Lumen	2.6 mm (.10 in)	2.6 mm (.10 in)
Length	758 mm (29.8 in)	758 mm (29.8 in)